
OBSOLETE **RETROWASH**

INSTRUCTIONS FOR USE

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BRUGSANVISNING

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BRUKSANVISNING

Any serious incident that occurs in relation to this device should be reported to the manufacturer and the competent authority of the Member State in which the user is established.

OBsolete

Instructions For Use

The Retrowash has been designed to facilitate intra-operative colonic lavage utilising retrograde irrigation. It is available with an outer diameter of (A) 35mm and (B) 46mm.

Components:

Retrowash Device - Flexible handpiece with three ports as follows:

- (1) Corrugated opening for attaching to the colon proximal to lesion.
- (2) Irrigation administration port sealed by a diaphragm with nipple.
- (3) Waste Outlet connected to a 1.6m clear flexible tube.
- (4) **Cubitainer Waste Receptacle** - Two detachable 10-litre containers with sealing caps for connecting to the Retrowash waste tubing.

- (5) **Giving Set** & (6) **Extension Set** - 2.0 metre Giving Set for connecting to irrigation solution, and 1.5m Extension Set for insertion through Retrowash and advancing along colon.

Preparation

1. Non-crushing intestinal clamps should be placed across the distal terminal ileum to prevent reflux into the small bowel. If necessary, manipulate the small bowel to deposit contents into the large bowel prior to clamping.
2. A non-crushing intestinal clamp should be placed across the colon proximal to the lesion to facilitate attachment of the Retrowash.
3. Select the appropriate Retrowash size. **The 46mm model should be used where the colon is noticeably distended, and the 35mm model should be used in all other cases.**

Set Up and Attachment

1. Ensure transected end of descending/sigmoid colon is at least 10cm distal to intestinal clamp.
2. Insert Retrowash handpiece until all corrugations are at least 2cm inside the colon. If the colon is too small and the product cannot be easily inserted, it should not be used.
3. Secure colon to the handpiece with the ties provided. The first should be applied around the handpiece **behind** the corrugations, then one between each of the corrugations.

4. Allow Retrowash waste tubing to fall to floor, and leave the Cubitainer (non-sterile) attached to end. **Do not discard the Cubitainer cap.**
5. Hang 3-litre bag of warmed irrigation solution of choice on drip stand. Close clamp on Giving Set and spike solution bag.
6. Connect Extension Set to end of Giving Set. **Apply sterile lubricant to bullet nose.**
7. Cut nipple of diaphragm **behind** marked line (i.e. distal to handpiece) and insert the bullet nose plus 1-2cm of the Extension Set. The Retrowash is now ready for use.

How to Use

1. Once satisfied that device and irrigation is securely attached, remove proximal intestinal clamp.
2. Irrigation should be in bursts of 250-500ml to avoid distension, with each burst followed by gentle manipulation of the colon. The free-flow irrigation rate is approximately 56 seconds/litre for 3-litre bags.
3. Release clamp on Giving Set and give first burst of irrigation. Gently manipulate the colon to assist breakdown of faecal material and enable drainage back through the Retrowash and down to the Cubitainer.
4. Advance Extension Set along the cleared colon (guide tubing externally by hand) to the next faecal compaction. Apply second burst of irrigation and manipulate as above.
5. Repeat steps 3 & 4 until the entire colon has been irrigated and the effluent runs clear through the waste tube. If necessary, suction may be connected to the Extension Set to remove residual effluent.
6. In the event of intussusception during irrigation, pressure may be released by disconnecting the Extension Set from the Giving Set.
7. Carefully withdraw Extension Set (now contaminated), and if desired insert a colonoscope through the same administration port and advance along the colon in the usual way.

NB: The remainder of the nipple should be cut off before inserting scope.

Removal & Disposal

1. Apply crushing clamps across colon proximal to connection with handpiece and divide.
2. Discard the tissue still attached to the handpiece.
3. Empty any residual irrigation fluid from waste pipe into Cubitainer, detach waste pipe, reseal Cubitainer with cap. Solidifying solution may be added.
4. Dispose of all parts as clinical waste. The maximum weight of the Cubitainer when full is 10kg, so 2 yellow clinical waste bags should be 'doubled up' when disposing.
5. Should the Cubitainer fill up before end of irrigation procedure, cease irrigating, detach and seal first Cubitainer as per point 3 above, and attach second Cubitainer to waste tubing.

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INTERMARK MEDICAL INNOVATIONS LTD